



**ALPHA2** **EQ**™

# Equine Alpha-2 Macroglobulin Safety Study and Clinical Observations

Brought to you by

**ASTARIA**™



## Equine Alpha-2 Macroglobulin Safety Study and Clinical Observations

Osteoarthritis (OA) and inflammatory joint disease are both mediated by numerous biomechanical and biochemical processes. Inflammatory cytokines have been demonstrated to increase production of metalloproteases that degrade cartilage, and the catabolic cartilage products further stimulate production of additional pro-inflammatory cytokines. Recently, increased attention has been given to inhibiting metalloproteases involved in cartilage catabolism, with the intention of modulating cartilage breakdown and preventing the cascade of inflammatory mediators involved in disease progression.

Alpha-2 Macroglobulin (A2M) has emerged as a unique potential treatment of cartilage-based pathology and inflammatory arthritides. The process of concentrating A2M for the treatment of a variety of joint and disc related issues in human medicine is well documented. A proprietary kit has been developed and FDA approved. Autologous concentrated A2M from plasma is currently in use to successfully treat various painful arthritides, including mild to moderate OA, PTOA, enthesopathies, and spinal discogenic back pain in the human market. A2M can not only inhibit the associated inflammatory cascade but also disrupt the catabolic process of cartilage degeneration.

This paper is to document the use of the Astaria Global proprietary kit (Alpha2EQ) system to concentrate equine A2M and utilize it in the treatment of a variety of orthopedic issues. The human kit is easily adopted to the equine patient and the process for preparing autologous A2M is identical.

Two important biochemical networks that contribute to OA pathology include pro-inflammatory cytokines and matrix metalloproteinases (MMPs).<sup>7</sup> A2M is a major plasma glycoprotein best known for its ability to inhibit a broad spectrum of serine, threonine and metalloproteases by a unique bait and trap method (Fig. 1).

A2M uses a 39–amino acid bait region that, when cleaved by a protease, induces a large irreversible conformational change that physically traps the protease within a steric cage. As part of the entrapment, the protease forms a covalent bond with A2M, exposing a receptor recognition site, triggering the endocytosis and eventual clearance of the A2M-protease complex. A2M has also been demonstrated to bind pro-inflammatory cytokines, such as TNF- $\alpha$  and IL-1 $\beta$ , and reduces the cytokine-induced up-regulation of collagenases in chondrocytes (Fig. 2).

Fig. 1. Pictorial representation of the mechanism by which A2M traps proteases. Each dimer traps a single protease. After the second trapped protease, the molecule then undergoes active transport for its elimination.

$\alpha_2$ -Macroglobulin 911





Fig. 2. A2M is chondroprotective against inflammatory cytokines. Treatment of bovine cartilage explants (blue column) with pro-inflammatory cytokines TNF- $\alpha$  and IL-1 $\beta$  (red column) induces chondrocytes within the cartilage to produce or activate proteases resulting in increased production of sulfated glycosaminoglycan (sGAG). Treatment with purified human A2M (green column) potentially inhibited cartilage catabolism. BCE, bovine cartilage explant.



Downloaded from ClinicalKey.com at UCLA -Cedars Sinai Medical Center November 05, 2016.  
For personal use only. No other uses without permission.  
Copyright ©2016. Elsevier Inc. All rights reserved.

## Study

For the purpose of establishing safety of the process in the horse, a pilot project utilizing 4 non-lame recipient mares was performed. Each mare was aseptically prepped for a blood draw. Blood was then processed using the kit provided. The medial femorotibial joint of each mare was injected with 15cc of the end product. Patients were evaluated for signs of inflammatory reaction and lameness at 2 hours, 12 hours, 24 hours and daily for the next 3 days. At no time was there any evidence of any reaction associated with the injection.

Based on the pilot project results, a clinical safety study was undertaken. A total of 50 clinical patients were chosen to receive joint injections for a variety of orthopedic conditions that would have traditionally been treated with existing modalities. The clinical study was performed as a multi-clinic study. Five equine hospitals<sup>1</sup> participated in the study and each contributed patients for evaluation of safety as well as clinical response. No controls were utilized, and the results are reported only as reaction to injection and positive or negative clinical response. A positive response was considered if a patient improved more than a full grade of lameness following treatment. In several instances, a second injection was required to achieve a positive clinical response. Patients that required a second injection were all treated with frozen product, which is indicated in the Clinical Response column. All patients were evaluated at 24 hours, 1 week and 2 weeks to determine adverse reactions and clinical response.

Signalment	Site Injected	Adverse Reaction	Clinical Response
<b>QH-Roping</b>	Carpus -Stifle	None	Positive-Received second injection at 6 months
<b>QH- Trail</b>	Coffin Joints	None	Positive
<b>QH-Cow Horse</b>	Coffin Joints-Stifle	None	Positive- Received second injection at 6 months
<b>QH-Cutting</b>	Sifles-MFT	None	Positive on 1
<b>QH-Roping</b>	Coffin Joint	None	Positive
<b>QH-Barrel Racing</b>	Carpus - TMT's	None	Positive
<b>Wmblood-Dressage</b>	Origin LF suspensory	None	Positive
<b>QH-Cow Horse</b>	Coffin Joints-Bursas	None	Positive-Received second injection at 1 month
<b>QH-Roping</b>	LH Coffin joint-Bursa	None	Positive

Signalment	Site Injected	Adverse Reaction	Clinical Response
<b>QH-Roping</b>	Coffin Joint	None	Positive
<b>QH-Barrel Racing</b>	Front Tendon Sheaths	None	Negative
<b>QH-Cow Horse</b>	Carpus	None	Positive
<b>QH-Roping</b>	Stifle-MFT	None	Required second injection
<b>QH-Barrel Racing</b>	Sacroiliac L	None	Positive
<b>QH-Ranch</b>	Coffin Joints-Bursas	None	Positive
<b>QH-Reining</b>	Coffin Joints	None	Positive
<b>Thoroughbred</b>	Fetlocks	None	Positive
<b>QH-Barrel Racing</b>	Coffin Joints	None	Positive
<b>Wmblood-Jumping</b>	Coffin joint-bursa	None	Positive- Received second injection at 6 months
<b>Arabian Show</b>	Coffin Joints	None	Pending
<b>QH-Roping</b>	Coffin joint-bursa	None	Required second injection
<b>QH-Roping</b>	Carpus	None	Positive
<b>QH-Barrel Racing</b>	RH suspensory	None	Positive
<b>QH-Barrel Racing</b>	Carpus	None	Positive
<b>QH-Roping</b>	Coffin joint-bursa	None	Required second injection
<b>QH-Barrel Racing</b>	Coffin Joints	None	Required second injection
<b>QH-Roping</b>	Coffin Joints-Stifles	None	Positive
<b>QH-Roping</b>	Navicular Bursa R&L	None	Positive- Received second injection at 8 months

Signalment	Site Injected	Adverse Reaction	Clinical Response
<b>Arabian Show</b>	Coffin Joints-LF Fetlock	None	Positive
<b>Warmblood</b>	Cervical facets C5-C7	None	Positive
<b>Arabian Hunter</b>	R&L Front Fetlocks-Coffin joints	None	Positive
<b>Warmblood</b>	LF Coffin Joint	None	Positive
<b>QH-Roping</b>	Navicular Bursa R&L	None	Positive
<b>Danish WB-dressage</b>	Coffin Joints	None	Positive
<b>Percheron X-pleasure</b>	RF coffin joint and bursa	None	Positive
<b>Belgian WB-eventing</b>	Coffin Joints	None	Positive
<b>ISH-Foxhunting</b>	Coffin Joints-TMT's	None	Positive
<b>Thoroughbred-eventing</b>	Coffin Joints	None	Positive
<b>Hanoverian X-Dressage</b>	Coffin Joints-TMT's	None	Positive
<b>Thoroughbred-Eventing</b>	Coffin Joints-Bursas - LF digital sheath	None	Positive
<b>Thoroughbred-Eventing</b>	Coffin Joints	None	Positive
<b>ISH-Eventing</b>	RF coffin joint	None	Positive
<b>Warmblood-Jumper</b>	RF coffin joint and bursa	None	Positive
<b>QH-Cutting</b>	TMT's	None	Positive
<b>Warmblood-Jumper</b>	Coffin joint-bursa RF	None	Positive
<b>Warmblood-Hunter</b>	L Inter-carpal joint	None	Pending

## **Conclusions:**

This study indicates that the application of Astaria Global proprietary kit (Alpha2EQ) system for concentration of A2M is safe and seemingly efficacious in the equine patient. No horse experienced any adverse reaction or joint flare post injection. In addition, the several patients that received frozen product injections did not experience any joint flares. The clinical response to frozen product injected was no different than that from freshly processed product.

Research is currently underway at Boise State University. Research to demonstrate the concentration levels of A2M, and the cytokine profiles present in the equine product has been completed and is under evaluation. Studies to determine the molecular size and configuration of the A2M molecule in the horse and other biologically active proteins that may be present is underway.

Idaho Equine Hospital - Dr. Stuart Shoemaker  
Virginia Equine Imaging - Dr. Kent Allen, Dr. Abigail Velting  
Brazos Valley Equine Hospital - Dr. Kurt Heite  
Kleinpeter Equine Hospital- Dr. Chat Kleinpeter  
Northwest Equine Performance- Dr. Mark Revenaugh